

AD-A189 380

DTIC FILE COPY

AD-.....

METABOLISM, MASS SPECTRAL ANALYSIS AND MODE  
OF ACTION OF TRICHOPELENE MYCOTOXINS

ANNUAL REPORT

15 JULY 1986 THROUGH 14 JULY 1987  
15 SEPTEMBER 1987

CHESTER J MIROCHA  
ROBERT J PAWLOSKY  
ROLAND GUNTHER

DTIC  
ELECTE  
JAN 26 1988  
S D

SUPPORTED BY

U.S. ARMY MEDICAL RESEARCH AND DEVELOPMENT COMMAND  
FORT DETRICK, FREDERICK, MARYLAND 21701-5012

CONTRACT NO. JAMD17-85-C-5204

UNIVERSITY OF MINNESOTA  
DEPARTMENT OF PLANT PATHOLOGY  
ST. PAUL, MN 55108

APPROVED FOR PUBLIC RELEASE; DISTRIBUTION UNLIMITED

The findings in this report are not to be construed as an official  
Department of the Army position unless so designated by other  
authorized documents.

20030128255

88 1 20 010

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER	2. GOVT ACCESSION NO.	3. RECIPIENT'S CATALOG NUMBER
4. TITLE (and Subtitle) Metabolism, Mass Spectral Analysis and Mode of Action of Trichothecene Mycotoxins		5. TYPE OF REPORT & PERIOD COVERED Annual Report 15 July 1986- 14 July 1987
		6. PERFORMING ORG. REPORT NUMBER
7. AUTHOR(s) Chester J. Mirocha, Roland Gunther and Robert J. Pawlosky		8. CONTRACT OR GRANT NUMBER(s) DAMD17-85-C-5204
9. PERFORMING ORGANIZATION NAME AND ADDRESS Department of Plant Pathology University of Minnesota St. Paul, MN 55108		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS 63763A 3M263763D807 AH.046
11. CONTROLLING OFFICE NAME AND ADDRESS USAMRDC FORT DETRICK FREDERICK, MD 21701-5012		12. REPORT DATE 15 September 1987
		13. NUMBER OF PAGES 59
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office)		15. SECURITY CLASS. (of this report) Unclassified
		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE
16. DISTRIBUTION STATEMENT (of this Report) Approved for public release; distribution unlimited.		
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)		
18. SUPPLEMENTARY NOTES		
19. KEY WORDS (Continue on reverse side if necessary and identify by block number)		
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) An analytical method for T-2, HT-2 and T-2-tetraol using multiple reaction monitoring (MRM) was initiated. The appropriate parent and daughter ions were selected. The sensitivity of the method is one part per billion. An isomer of HT-2 was detected. The mass spectra of H-1 (parent and daughter ions) were obtained. Toxic Fusarium isolates were identified.		

## SUMMARY

The development of an analytical method for the analysis of T-2, HT-2 and T-2-tetraol in blood and urine was initiated. The objective is to use multiple reaction monitoring (MRM) facility of the VG 7070EQ tandem mass spectrometer. The advantage of the latter is that the noise from the biological matrix is kept at a minimum and as a result the sensitivity is increased. Essentially, the magnet is parked on the parent ion and only this fragment is allowed to pass through the first mass spectrometer. After the selected fragment passes into the quadrupole #1, it is decomposed into fragments by collision activated (in argon plasma and 20 30 volts) decomposition. The daughter ions are detected by quadrupole #2. Using this method, T-2 and HT-2 can be detected at one part per billion. *Keywords: T-2 Toxin, HT-2 Toxin*

In the development of the MRM method of analysis for T-2-TFA, we discovered an interfering parent fragment of 478 that had a retention time close to that of T-2. The substance did not present a problem because its daughters were different from that of T-2-TFA. This also illustrates that although the parents of a substance are the same, the daughters are not.

The fragmentation map of the TFA derivatives of T-2, HT-2 and T-2-tetraol as well as their daughter ion spectra have been obtained.

An isomer of HT-2 (C-4 acetate) was discovered in a Fusarium culture using the daughter ions of HT-2 as a probe. Its structure was determined using parent-daughter mass spectra. The detection of the isomer in the Fusarium culture, composed of a vast mixture of compounds, illustrates the utility of the daughter ion library as a practical analytical probe.

The daughter ion fragments for T-2, HT-2 and T-2-tetraol were determined using respectively parents  $m/z$  478, 532, and 568. T-2-tetraol can also be detected by using the parent 330 which yields the daughter 216.

This report also lists all the Fusarium species isolated from Norway and New Zealand and the toxins they produce. Those from Norway almost exclusively produce the hemorrhagic toxin called H-1. This toxin was also found in New Zealand on the South Island.



Accession For	
NTIS CRA&I	<input checked="" type="checkbox"/>
DTIC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By	
Distribution /	
Availability Codes	
Dist	Avail and/or Special
A-1	

## FOREWORD

Citations of commercial organizations and trade names in this report do not constitute an official Department of the Army endorsement or approval of the products or services of these organizations.

# TABLE OF CONTENTS

	PAGE
Summary . . . . .	i
Foreword. . . . .	ii
Multiple Reaction Monitoring (MRM) of T-2 and HT-2 toxins and T-2-tetraol. . . . .	1
Analysis of HT-2 toxin. . . . .	1
Detection of an isomer of HT-2 toxin (TFA) using the daughters of 180, 138 and 121 . . . . .	1
Analysis of T-2-tetraol . . . . .	2
Conditions used in MRM analysis . . . . .	2
Structure Elucidation of the C-4 acetate isomer of HT-2 . . . . .	2
FIGURES:	
Figure 1. Mass spectrum (E.I.) of the T-2 TFA derivative and fragmentation map of HT-2 toxin. . . . .	4
Figure 2. Daughter ion spectrum of HT-2 (TFA) using the M+ 616. . . . .	4
Figure 3. Daughter ion spectrum of m/z* 532 of HT-2 (TFA) . . . . .	5
Figure 4. Daughter ion spectrum of m/z+ 514 of HT-2 (TFA) . . . . .	5
Figure 5. Daughter ion spectrum of m/z+ 455 of HT-2 (TFA) . . . . .	6
Figure 6. MRM analysis of urine spiked with HT-2 at 50 ppb using parent 532 and daughters 180, 138 and 121. . . . .	6
Figure 7. Mass spectrum of T-2-tetraol (TFA) and fragmentation map of daughters from M+682. . . . .	7
Figure 8. Daughter ion mass spectrum of 682 of T-2-tetraol (TFA). . . . .	7
Figure 9. Daughter ion mass spectrum of 569 of T-2-tetraol (TFA). . . . .	8
Figure 10. Daughter ion mass spectrum of 568 of T-2-tetraol (TFA) . . . . .	8
Figure 11. Daughter ion mass spectrum of 455 of T-2-tetraol (TFA) . . . . .	9
Figure 12. Daughter ion mass spectrum of 330 of T-2-tetraol (TFA) . . . . .	9
Figure 13. MRM analysis of urine spiked with T-2-tetraol (TFA) at 50 ppb using 330 and 568 as parents . . . . .	10
Figure 14. Structure of fragments of the C-4 isomer of HT-2 (TFA) toxin . . . . .	11
Figure 15. Continuation of the structures started in Figure 14. . . . .	11
Figure 16. Positive chemical ionization mass spectrum in methane of the C-4 isomer of HT-2 (TFA) toxin . . . . .	12
Figure 17. Daughter ion spectrum of the C-4 isomer of HT-2 (TFA) using 532 as the parent . . . . .	12
Figure 18. Daughter ions of C-4 acetate isomer of HT-2, C-4 isomer (acetate) of T-2 and C-3 isomer (acetate) of T-2 . . . . .	13
Figure 19. Positive chemical ionization mass spectrum of T-2 (TFA) and HT-2 C-4 acetate isomer of HT-2 (TFA) . . . . .	13
Implementation of T-2 Toxin Analysis Using Multiple Reaction Monitoring . . . . .	14
Structure Identity of H-1 Toxin . . . . .	14
Figure 1. Electron impact mass spectrum of H-1 toxin . . . . .	15
Figure 2. Chemical ionization mass spectrum of H-1 in methane. . . . .	15
Figure 3. Ultra violet absorbtion spectrum of H-1. . . . .	16
Figure 4. Infrared absorbtion spectrum of H-1. . . . .	17

# TABLE OF CONTENTS (Continued)

	PAGE
Table 1. Results of combustion analysis of H-1 by Galbraith. . .	18
Pathologic Studies of H-1 toxin on rats . . . . .	19
Table 1. LC50 values of H-1 toxin on rats, mice and guinea pigs. . .	22
Table 2. Weight gain and feed consumption of rats given H-1 in their diet at 4 ppm. . . . .	23
Mass Spectral Library of Trichothecenes and Other Mycotoxins	
Produced by <u>Fusarium</u> . . . . .	24
Mass Spectral Library of Trichothecenes. . . . .	24
Use of Library for Identification Purposes . . . . .	25
Figure 1. Resolution of crude extract of TDP-1 producing isolate of <u>Fusarium</u> . . . . .	26
Figure 2. Report of library search on culture shown in Fig. 1. .	26
Figure 3. Comparison of unknown in <u>Fusarium</u> culture with library entry . . . . .	27
Figure 4. Mass spectral library search of scan #245 suspected of being TDP-2 (acetate of TDP-1). . . . .	27
Figure 5. Comparison of scan #245 with library entry . . . . .	28
Figure 6. Comparison of the total ion chromatogram of the <u>Fusarium</u> extract with m/z+ 218 (base peak of TDP-1) .	28
Interfering Substance in Human Blood and its Analysis by the Daughter Ion Mass Spectral Library . . . . .	29
Figure 7. Daughter ion mass spectra of parent m/z+478 of the interfering substance found in blood. . . . .	30
Figure 7-A. Resolution of the interfering substance found in blood from T-2 toxin (TFA). . . . .	31
Figure 8. Results of library search of the interfering substance with the daughter ion library . . . . .	32
Figure 9. Results of the daughter ion mass spectral library for T-2 toxin . . . . .	32
Figure 10. Comparison of the daughter ion library spectrum of T-2 (TFA) with the suspected T-2 found in the urine sample . . . . .	33
Figure 11. Parameters used in the library search routine. . . .	33
Analysis of T-2 (TFA) Using Multiple Reaction Monitoring. . . . .	34
Figure 12. Full scan mass spectrum of T-2 (TFA) in E.I. using a source temperature of 150 C. . . . .	35
Figure 13. Daughter ions of m/z+ 401 of T-2 (TFA). . . . .	35
Figure 14. Daughter ions of m/z+ 478 of T-2 (TFA). . . . .	36
Figure 15. Analysis of 2 ppb T-2 in blood by MRM. The daughter at m/z+ 180 is shown . . . . .	36
Figure 16. Analysis of 2 ppb of T-2 in blood by MRM. Shown are daughters 180, 138, 121 and the parent at 478. . . .	37
Daughter ions of the TFA Derivatives of T-2, HT-2 and T-2-tetraol	38
Figure 17. Electron impact mass spectrum of T-2 (TFA) and its fragmentation map of daughter ions . . . . .	39

## TABLE OF CONTENTS (Continued)

	PAGE
Figure 18. Electron impact mass spectrum of HT-2 (TFA) and its fragmentation map of daughter ions . . . . .	39
Figure 19. Electron impact mass spectrum of T-2-tetraol (TFA) and its fragmentation map of daughter ions . . . . .	39
 Analysis of T-2 Toxin in Human Blood Using Tandem Mass Spectrometry-Gas Chromatography . . . . .	
Amending Blood with T-2 and Extraction. . . . .	40
Derivatization . . . . .	40
GC Conditions. . . . .	40
Source Conditions. . . . .	40
MS/MS Tuning . . . . .	41
Collision Activated Decomposition. . . . .	41
Multiple Reaction Monitoring and Results . . . . .	41
Figure 1. On-Column injection device . . . . .	44
Figure 2. On-Column injection of T-2 (TFA) and E. I. spectrum. . . . .	45
Figure 3. Daughter ions obtained from parent 478 of T-2 (TFA). . . . .	45
Figure 4. Parameters used in MRM analysis. . . . .	46
Figure 5. Ratios and retention time of daughters 180, 138 and 121 of T-2 (TFA) analyzed by MRM. . . . .	46
Figure 6. On-Column injection of an extract of blood amended with 2 ppb T-2 (TFA) toxin. . . . .	47
Figure 7. Daughter of parent 478 found in the extract of blood. The latter's retention time was different from that of T-2 (TFA). . . . .	47
Figure 8. Resolution of 2 ppb T-2 toxin in a blood matrix by splitless capillary column introduction as compared with on-column injection. . . . .	48
Figure 9. Resolution of the 180 daughter of T-2 (TFA) in a blood matrix using the splitless method of introduction . . . . .	48
Figure 10. Resolution of daughter 121 of T-2 (TFA) at 2 ppb in blood using the splitless method of introduction . . . . .	49
Figure 11. Resolution of Daughter 138 of T-2 (TFA) at 2 ppb in blood using the splitless method of introduction . . . . .	49
Figure 12. T-2 (TFA) standard analyzed by MRM (splitless). . . . .	50
 Toxicity and Toxins Produced by Isolates of <u>Fusarium</u> Species Obtained from Norway and New Zealand . . . . .	
Table 1. Summary table of toxicity of <u>Fusarium</u> isolates to rats. . . . .	51
Table 2. Production of mycotoxins of <u>Fusarium</u> collected in Norway and New Zealand after growth on rice. . . . .	53
Table 3. Summary of mycotoxins detected in extracts of <u>Fusarium</u> used in the toxicity studies . . . . .	54
Table 4. Mycotoxins found in fifteen toxic isolates of <u>Fusarium</u> collected in New Zealand . . . . .	56
Distribution List . . . . .	57
	59

MULTIPLE REACTION MONITORING (MRM) WITH T-2 TOXIN, HT-2 TOXIN  
AND T-2-TETRAOL

In Quarterly Report number VI, we reported on the analysis of T-2 toxin in a blood matrix using the MRM procedure with the VG7070EQ instrument. We were able to obtain what should be the maximum in selectivity and sensitivity of the scanning portions of the instrument. MRM is a tandem mass spectral scanning technique that forces the electro static analyzer and analyzers of the magnet and quadrupole to focus on specific reactions (parents to daughter ions) and thereby increase sensitivity and selectivity. We report here the expansion of this technique to include HT-2 toxin and T-2-tetraol.

Selectivity of MRM relies on specific parent and daughter ion relationships. Parent ions of a certain molecule will produce unique unambiguous daughter ions of precise abundance using a specified set of reaction conditions. The daughter and parent ions taken together make up what is known as a fragmentation map. The latter is a pathway of decomposition of the ionized molecule and suggests the structure of the molecule. They are also valuable analytically. The majority of the total ion current will reside within only a few of the fragments usually found in the low mass region. Normally the low mass region would have many interfering substances from the biological matrix but in the case of MRM, this fault is eliminated. The reason for this is that only a single parent is allowed through the analyzer and on to the collision chamber where the daughter ions are formed. Thus the MRM procedure is well suited to take full advantage of high abundance ions as well as highly specific reactions.

ANALYSIS OF HT-2 TOXIN: The fragmentation map and daughter ion spectra of HT-2-(TFA) is presented in figures 1-5. Most of the abundance of the ion current is located in fragments 121, 138 and 180. These ions are daughters of parent 532. They are presumed to have the same structure as the same daughters found in T-2-TFA but derived from parent ion 478. Both parent ions ( 532 for HT-2 and 478 for T-2) are formed from a loss of 84 mass units from the molecular ion. Since these are high abundance ions and are formed from a significant precursor of the molecules intact skeleton, they are an excellent choice for MRM analysis. Figure 6 shows the results of analysis of HT-2 from urine using daughters 180, 138 and 121. Note the prominent ions and lack of noise from the biological matrix.

DETECTION OF AN ISOMER OF HT-2 TOXIN(TFA) USING THE DAUGHTERS 180, 138 AND 121: The fragments of an isomer of HT-2 TFA is shown in figure 14 and 15. The isomerism resides around C-4 carbon and



was detected as shown in the caption of figure 14. The diagnostic ion for the C-4 isomer of HT-2 toxin is  $m/z+227$ . Figure 16 shows the positive chemical ionization mass spectrum of iso-HT-2. The base peak of the most abundant HT-2 isomer (C-15 acetate) is 455. However, the base peak of the iso-HT-2 is 401 just as found in T-2 toxin (figure 19). The fragment 401 is formed by loss of the isovaleroxy group (101) to yield 515 followed by a loss of the C-15 TFA group (114) to yield 401.

ANALYSIS OF T-2-TETRAOL: The daughter ion fragmentation map of T-2-tetraol (figure 7) is different from that described for T-2 and HT-2 toxins. They differ in ion abundance as well as the structure of the primary and secondary fragments. The two most significant reactions are the parent ion 569 going to daughter 455 and parent 330 going to 216. The latter are abundant, reproducible, give the minimum of background signal and are easy to monitor. The MRM analysis of T-2-tetraol extracted from urine is shown in figure 13. The parent ion 330 yields a very prominent daughter in 216. On the other hand, parent 568 yields a very poor daughter at 454. The parent of choice is  $m/z+569$ .

CONDITIONS USED IN MRM ANALYSIS: The MRM analyses were done using two reaction groups. The first one covered the retention time of T-2-tetraol and the second covered T-2 and HT-2. The dwell time was 100 milliseconds, with an analyzer switch time of 20 milliseconds. The mass spectrometer tuning and calibration procedures were as described by the manufacturer i.e. calibrate the front end first followed by the quad. The electron voltage was 70, collision cell voltage was 14-16 volts, the collision cell argon pressure was 10<sup>-6</sup> torr.

The derivatives used were trifluoroacetate made from trifluoroacetic acid anhydride. The capillary column was a DB-5 (J&W Scientific) 0.25 mm id. The column temperature and program was 80 to 300 degrees C. at 25 degrees per minute.

The urine was spiked with an equivalent to give 50 parts per billion. Samples were extracted according to procedures reported in earlier reports using BONDELUTE columns.

#### STRUCTURE ELUCIDATION OF THE C-4 ACETATE ISOMER OF HT-2 TOXIN:

The C-4 isomer is called iso-HT-2 analogous to isoT-2 toxin. It is 3,15-dihydroxy-4-acetoxy-8-[3-methyl-butyloxy]-12,13-epoxy-trichothec-9-ene. The structure was determined by deduction using the mass spectrum and daughter ions formed. The isomer was found in a culture of a *Fusarium* isolate grown on rice in the laboratory. The primary mass spectra of the TFA derivative in both electron impact and chemical ionization was confirmed as

being 616 (figure 16). The mass spectra were similar but not identical to the C-15 isomer of HT-2 and the C-4 isomer of T-2 toxin. The daughter ion spectra are helpful in identifying the isomer because of the sequence of loss of the isovaleroxy group at C-8 (loss of 101) followed by the loss of the C-15 acetate as in the C-15 isomer (loss of 60) or loss of 114 amu in the case of the C-4 isomer. Thus the C-4 isomer should yield a either  $m/e+515$  or  $532$  followed by either  $401$  or  $413$  and finally resulting in  $m/z+227$ . Iso-T-2 toxin, in contradistinction, will yield an  $m/z+287$ .

The C-15 acetate isomer of HT-2 toxin should yield a  $m/z+455$  because of the loss of an acetate (60amu) from  $m/z+515$ . The latter is formed after loss of 101 amu from the isovaleroxy group on C-8 (figure 19). The C-4 acetate isomer will lose 114 amu units from the TFA formed at C-15 and yield  $m/z+401$  as in T-2 toxin.

FIGURE 1. Mass spectrum of T-2 TFA derivative. The  $M^+$  is 616 and the origin of daughters is as shown in the fragmentation scheme. Thus  $m/z$  616 gives rise directly to 532 which in turn is the parent of  $m/z$  fragments 180, 138 and 121. These are the same daughters as those found in the parent fragment  $m/z$  478 of T-2-TFA. The latter daughter fragments are very intense and are the choice for both quantitative and qualitative determination of HT-2-TFA. This is an example where two closely related derivatives have the same daughters but different parents i.e. 532 in HT-2 and 478 in T-2.




FIGURE 2. Daughter ion mass spectrum of HT-2-TFA using the  $M^+$  616 as the parent fragment. Note that  $m/z$  180, 138 and 121 are present in the spectrum but vary in intensity from that of the full scan E.I. mass spectrum shown in figure 1. Daughters 556, 423, 333, 264 and 220 are not found in the full scan mass spectrum indicating that different reactions are occurring under E.I. conditions. Note that  $m/z$  455, a classic indicator of HT-2 and its derivatives is not found as a daughter of 616.




FIGURE 3. Daughter ion mass spectrum of parent  $m/z$  532 of HT-2-TFA using argon as the collision plasma. Daughters 121, 138 and 180 are the products of collision activated decomposition of 532. The latter are very intense and stable, hence useful in MRM experiments.




FIGURE 4. Daughter ion mass spectrum of HT-2-TFA using  $m/z$  514 as the parent. Note that fragment 454 occurs as the daughter which is analogous to  $m/z$  455 of the full scan E.I. spectrum. It is presumed that 455 can be formed just as well as 454 but it depends on which reaction species the hydrogen follows.




FIGURE 5. Daughter ion mass spectrum of parent 455 of HT-2-TFA. Note that  $m/z+455$  is a stable species with little activity in terms of daughter fragments. This suggests that this fragment has an unsaturated ring structure which normally confers stability i.e. double bonds in the ring are stable.





FIGURE 6. Analysis of HT-2-TFA toxin extracted from urine, spiked at 50 ppb, using multiple reaction monitoring. Parent fragment 532 was allowed to pass through the magnet portion of the MS after which it was bombarded with electrons in the presence of argon in the first quad region and scanned for daughters in the second quad region. Boxes B thru D represent the total ion current of  $m/z+$  fragments 180, 138, and 121 respectively.






FIGURE 7. Electron impact mass spectrum of the TFA derivative of T-2-tetraol. The scheme of daughter ions and their parents is presented in the right hand corner. The major diagnostic fragments are  $m/z$  568, 455, 330 and 217.

FIGURE 8. Daughter ion mass spectrum of fragment 682 ( $M^+$ ) of the TFA derivative of T-2-tetraol. The most important daughter fragment ion is  $m/z$  568 which is the base peak of the E.I. spectrum




FIGURE 9. Daughter ion mass spectrum of  $m/z+569$  giving rise to  $m/z+455$  and 328.




FIGURE 10. Daughter ion mass spectrum of parent fragment 568 and the resulting daughter fragments. Note that  $m/z+455$  (daughter of 569) is not present in this spectrum.




FIGURE 11. Daughter ion mass spectrum of the TFA derivative of T-2-tetraol using  $m/z$  455 as the parent. Although there are a number of low intensity daughters, none are of interest in the use of MRM.




FIGURE 12. Daughter ion mass spectrum of the TFA derivative of T-2-tetraol parent fragment 330. The resulting daughter is 216 which is an important and useful fragment for MRM. We recommend its use for MRM of T-2-tetraol.






FIGURE 13. Analysis of human urine spiked with T-2-tetraol (TFA) at 50ppb and analyzed by MRM using both 454 and 330 as parent ions. Box A shows the total ion current obtained when  $m/z$  568 is used as the parent. Note the noisy base line which precludes the use of this fragment. On the other hand, box B shows the results of 330 as the parent and 216 as the daughter. Note the excellent base line with a signal to noise ratio of about 100/1 which suggests that the use of 568 as a parent in a biological matrix will give us sensitivity down to 2ppb.

UR48LA 8-APR-87 28:24 78ED Acft:RJP Sus:VERRA  
GR 1 A: 568.0000 B: 330.0000  
Text: SPIKED URINE MATRIX FOR T-2/MT2/T248L 50 PPB




FIGURE 14. Fragments of the TFA derivative of the C-4 isomer (acetate at C-4 instead of C-15) of HT-2 toxin. All of the above fragments and structures have been verified by daughter ion mass spectrometry. TFA-HT-2 toxin will yield a  $m/z+227$  derived from 616 minus the isovaleroxy group to yield  $m/z+515$  followed by loss of the TFA from C-15 yielding 401 and then a loss of 60 amu from the C-4 acetate to yield 341 and finally 114 amu from the loss of the C-3 TFA to give 227. The C-15 isomer will also yield  $m/z+227$  but from another parent ion. Finally, the C-3 isomer will yield  $m/z+287$  by loss of the C-8 group followed by loss of 114 amu each from C-15 and C-4 respectively. Fragments 227 and 287 are useful in detecting isomers of T-2, HT-2, TC-1 and TC-3.




FIGURE 15. Continuation of the breakdown scheme of HT-2 as shown in figure 14. Note that fragments 121, 138 and 160 are unsaturated rings and hence stable in the mass spectrometer. They are further degraded to  $m/z+166$  and 121.




FIGURE 16. Positive chemical ionization mass spectrum of iso-HT-2 toxin (acetate in the C-4 position). Note that the base peak is  $m/z+401$  instead of the traditional  $m/z+435$ . This is due to the fragmentation of the  $M^+$  at 516 which loses the C-8 isovalerolxy to create 515 and then the C-15 TFA ( $-114$ ) to form  $m/z+401$ . The latter mass is characteristic of TFA-T-2 toxin in PCI.




FIGURE 17. Daughter ion mass spectrum of the TFA derivative of HT-2 toxin (C-4 isomer) using  $m/z+532$  as the parent. This is compared with the daughters of  $m/z+532$  from TFA-T-2 toxin. Note the difference in daughters due mainly to the substitution of the TFA group on carbons 15, 4 and 3. The C-15 isomer loses 84 amu to form 532 and then 298 (the A ring and the epoxide plus C-3 and C-4) to form the unsaturated B ring at  $m/z+234$ . See figure 15




FIGURE 18. Daughter ion mass spectra (TFA derivatives) of the C-4 isomer of HT-2, C-4 isomer of T-2 and the C-3-T-2 toxin. All daughters have been generated from the same parent ( $m/z$  401). Although this is true, their structures are not the same. As an example, iso-T-2-TFA should yield  $m/z$  287 (middle box) and does whereas TFA-T-2 should yield  $m/z$  227 (bottom box) and does. The daughter ion method is excellent for differentiating this group.




FIGURE 19. Positive chemical ionization (methane) mass spectrum of TFA-T-2 toxin and TFA-HT-2(C-15 isomer). This figure shows that ion 341 of T-2 is formed from 563 by loss of an acetate (60amu) whereas in HT-2, ion 341 is formed by loss of a TFA group (114) from 616.

## IMPLEMENTATION OF T-2 TOXIN ANALYSIS USING MULTIPLE REACTION MONITORING.

The objective of this study is to develop a fool-proof method of analysis of human blood and urine samples for traces of T-2 toxin, HT-2 toxin and T-2-tetraol using multiple mass spectrometry. The reason for this study is to have on record a method to deal with chemical warfare agents such as the trichothecenes so that an unequivocal answer can be given as to their presence in biological samples.

Our criteria are to achieve: 1. ultimate sensitivity (1ppb) and 2. uncompromised detection by the elimination of biological matrix effects. We hope to surpass the results obtained by use of the Townsend Oxygen Discharge Electrode.

To date we have had good results down to 5ppb using a single daughter ion as the parent for detection of the characteristic fragments. We are studying the daughters in order to determine which of the many daughters lend themselves for good quantitation. A more complete report will be forthcoming in the next quarterly report.

## STRUCTURE IDENTITY OF H-1 TOXIN

The mass spectra in both E.I. and C.I. are shown in figures 1 and 2. The molecular ion in EI is 428 and this is confirmed by the molecular ion (428 plus 1) obtained in positive CI (methane). This was also the molecular ion obtained in analysis by Fast Atom Bombardment. Both the EI and CI spectra are rich in large fragment ions in the high mass region of the spectrum. Thus,  $m/z$  428, 385, 368, 354, 323 and 296 are very intense and informative ions useful in both interpretation of structure and in quantitation.

The ultraviolet absorption spectrum shows 210 and 250 millimicrons as absorption maxima. The UV spectrum does not appear to be too useful.

The infrared spectrum has intense peaks which will prove to be useful in interpretation of the structure after more detailed information is available from NMR analyses. Notable are the absorption maxima at 1750, 1650, 1200 reciprocal cm.

Combustion analysis revealed a ratio of carbon, hydrogen and oxygen of 64, 6 and 25% respectively.

Analysis by high resolution mass spectrometry suggest an empirical formula of  $C_{23}H_{24}O_8$ .

Structure analysis is being done by NMR. Large sample size preparations are being made in order to facilitate structure identification.

H1JUL290036 x1 Bgd=30 29-JUL-86 14:58-0:04:33 70E0 EI+  
 Bp=0 I=197mv Ha=0 TIC=15144000 SU Ront: Sys:RACAL  
 Text: H1 DIP 20 EV Cal: P1 #36 1.0  
 1291000  
 45000




Figure 1. Electron impact.

H1JUL290023 x1 Bgd=17 29-JUL-86 15:55-0:03:02 70E0 CI+  
 Bp=0 I=430mv Ha=0 TIC=45800000 Ront: Sys:RACAL  
 Text: H1 DIP METHANE CI 200 EV Cal: P1 #23 3.0  
 2800000




Figure 2. Chemical ionization.




Fig. 3. The UV spectrum of hemorrhagic factor H-1 in purified form (A) and rice inoculated with *F. oxysporum* N17B extracts (B) in ethanol.




Figure 4. Infrared spectrum of H-1 toxin (KBR pellet).



HARRY W. GALBRAITH, Ph.D.  
CHAIRMAN OF THE BOARD

KENNETH S. WOODS  
PRESIDENT

GAIL R. HUTCHENS  
EXECUTIVE VICE-PRESIDENT

VELMA M. RUSSELL  
SECRETARY TREASURER

# GALBRAITH

*Laboratories, Inc.*

QUANTITATIVE MICROANALYSES

ORGANIC - INORGANIC

KNOXVILLE, TENNESSEE 37921

P. O. BOX 4137  
2323 SYCAMORE DR.

PHONE 546-1335  
AREA CODE 615

Mr. Hamed K. Abbas  
University of Minnesota  
Plant Pathology / 495 Boylston  
1991 Buford Circle  
St. Paul, Minnesota 55108

August 25, 1986

Received: August 18th

Dear Mr. Abbas:

Analysis of your compound gave the following results:

Your #,	Our #,	% C,	% H,	% N,	% O,
H-1	P-9950	63.80	6.25	0.11	25.19

The molecular weight results will follow later.

This is confirming our telephone call of August 25, 1986.

Sincerely yours,  
GALBRAITH LABORATORIES, INC.

  
Gail R. Hutchens  
Exec. Vice-President

GRH:pd

## **Pathogenic Studies of H-1 Toxin on Rats**

We conducted acute trials with recrystallized H-1 toxin and a subacute trial with H-1 contaminated feed on rats. Acute trials with ethylacetate extracts of contaminated feed were performed on mice and guinea pigs as well.

### **Trial I: Gavage of Rats with Purified H-2 Toxin.**

As a pilot trial for studies with purified toxin, groups of 2 three week, female Sprague Dawley rats were gavaged with a single dose of 50, 25 or 4 mg/kg body weight recrystallized toxin in 0.5 ml 40% ethanol. Both animals in the 50 mg and 25 mg groups died within 24 hours, as did one of the animals at the 4 mg level. The other animal in the 4 mg group remained clinically normal and was euthanized at 24 hours. Gross lesions included severe hemorrhage of the mucosa of the glandular stomach in one animal of the 50 mg and 4 mg groups, severe cardiac hemorrhage in the rat of the 4 mg group and hyperemia of the duodenum in one rat of the 25 mg group and the rat of the 4 mg group. Histologic evaluation is still in progress, but preliminary results show that purified toxin causes similar lesions to those produced by crude toxin.

### **Trial II: Gavage of Rats with a Low Dose of Purified Toxin**

In an attempt to determine an effective dosage for a future acute study, 9 rats were gavaged with 4 mg/kg recrystallized toxin. Five of 9 died within 24 hours; the remaining animals remained clinically normal and were euthanized at 24 hours. There were no gross lesions.

### **Trial III: Subacute Feeding Trial in Rats**

In order to assess the effect of low level exposure to H-1 toxin, eleven 3 week old, female rats were fed a mixture of 10% contaminated

rice (4 ppm toxin) and 90% complete feed. Food and water intake and body weights were recorded at 48 hour intervals. Packed cell volumes, total serum protein levels, total and differential leukocyte counts were performed on 2 different animals at 48 hour intervals. Bone marrow smears, serum chemistry panels and necropsies were performed at 21 days. Ten controls were fed complete rat diet, ad libitum, for 15 days and 10 controls were fed approximately the same amount of food as was ingested by the treated animals. The trial was ended at 21 days, because several animals died at 20 days.

All treated animals became depressed and diarrheic, but all controls remained clinically healthy. Weight gain and feed conversion data are presented in Table 1. Complete results of serum chemistry, hematology, bone marrow evaluation and histologic examination are still pending, but the treated animals developed greatly elevated blood leukocyte levels.

**Table VI:** Administration of Crude Extracts to Mice and Guinea Pigs

As a pilot trial to evaluate the comparative toxicity of H-1, ethylacetate extracts of contaminated rice were administered by gavage to 180 g female guinea pigs and by intraperitoneal injection to 3 week, female mice. Groups of 2 guinea pigs were given a single dose of extracts of 10, 20 or 50 g of contaminated rice (44 ppm toxin) in 1-2 ml 20% ethanol. One animal of the 50 g group died within 24 hours with locally extensive gastric hemorrhage. The other animals were euthanized at 24 hours and did not develop clinical or morphologic abnormalities. Histologic examination is pending.

Groups of 3 three week, female mice were dosed with extracts of 0.63 and 1.25 g of contaminated rice. Two of the animals of the 1.25 g

group died within 24 hours, with one showing myocardial hemorrhage and both showing hyperemia of the duodenum. No clinical or gross abnormalities were found in the other mice; histologic examination is still pending.

#### Comments

These preliminary results indicate that acute administration of purified H-1 toxin results in similar clinico-pathologic effects to those produced by crude toxin. The lack of availability of purified toxin has hampered this portion of the study to some degree, but current production should resolve this problem.

Preliminary results also indicate that H-1 toxin produces lesions in mice and guinea pigs which are similar to those in rats. It also appears that guinea pigs are more resistant to orally administered toxin than rats.

The subacute feeding trial showed high toxicity at a level of approximately 4 ppm. Several animals died on day 20 and the others were very sick when euthanized. We have no explanation for the elevated blood leukocyte counts, but they may be a sequella to the lymphocidal effect noted with the pure toxin.

Our immediate plans are to complete the described studies and to perform more acute studies with purified material. In addition, the low level feeding trial will be repeated and the effects on the humoral response to foreign antigen (sheep red cells) and the in vitro lymphocytic response to mitogens will be assessed.

TABLE 1. LC50 VALUES OF THE HEMORRHAGIC FACTOR NAMED H-1 WHEN ADMINISTERED TO RATS, MICE AND GUINEA PIGS.

SPECIES	PURITY OF H-1	AGE DAYS	METHOD OF ADMINISTRATION	LC50 MG/KG	TOXIC SIGNS
RAT	crystals	20	ORAL	4	HI, HS, HT, HH, D
RAT	extract	20	ORAL	3	HB, HI, HS, HT, HH, D
MOUSE	extract	20	I.P.	4	HI, HH, D
GUINEA PIG	extract	20	ORAL	12	HS, D

All animals were purchased from BIO-LAB Corp. St. Paul, MN. Each was intubated or injected with toxin based on the gram equivalents of H-1 found in the toxic culture. Control animals received the same volume of 10, 20 or 40% aqueous ethanol as given to the experimental animals.

TOXICITY LEGEND:

HI= intestinal hemorrhage.  
 HS= hemorrhage in the stomach.  
 HT= hemorrhage in the thymus.  
 HH= hemorrhage in the heart.  
 D= death.

**Table 2. Weight gain and feed consumption of rats given a chronic exposure to H-1 toxin in their diet.\***

	Treated		Controls		Pair-fed Controls	
	Day 15	Day 21	Day 15	Day 21	Day 15	Day 21
Average						
Wt Gain	-18.4	-14.0	57.8	--	56.0	67.4
(g)						
Wt Gain/	-0.2	--	0.5	--	0.5	--
Feed Consumption						

\*Rats were given a diet containing 90% nutritionally balanced feed and 10% Fusarium culture containing 4 ppm H-1.

MASS SPECTRAL LIBRARY OF TRICOTHECENES AND OTHER MYCOTOXINS PRODUCED  
BY FUSARIUM

I. One of the objectives of this contract is to construct a mass spectral library of trichothecenes so that it could be used in the identification of this group of toxins and also facilitate the study of metabolism of T-2 and other members of this group. A library consisting of approximately 250 mass spectra of trichothecenes was constructed on the Hewlett Packard 5987 Real Time Executer System and delivered to USAMRIID a few years ago. The latter was developed on the Hewlett Packard Quadrupole Mass Spectrometer. It is our intention to transfer these spectra to the VG 7070 system which appears to be a more versatile and efficient system. Moreover, we are also developing a library of the daughter fragments of the trichothecenes on the VG 7070 EQ. The latter is designed to give unequivocal identification of the trichothecenes. We are including a new dimension in this library in so far that it will contain other potent toxins produced by species of Fusarium, i.e. other than the trichothecenes.

IN SUMMARY WE ARE DEVELOPING:

MASS SPECTRAL LIBRARY OF FULL SPECTRA (VG SYSTEM)

MASS SPECTRAL LIBRARY OF DAUGHTER IONS (VG SYSTEM)

A. MASS SPECTRAL LIBRARY OF TRICOTHECENES (FULL SPECTRA)

An example of our current effort in constructing this library is enclosed. Of more importance, we are showing its utility in identification of a toxic substance produced by Fusarium and present in a biological matrix. An example of a search follows.

This library is called FUSTOX and presently contains the following spectra:

T-2 toxin

3-hydroxy-4,15,-diacetoxy-8-[3-methyl-butyryloxy]-12,13-epoxy-trichothec-9-ene (TFA)

T2-tetraol (3,4,15-tetrahydroxy-12,13-epoxytrichothec-9-ene(TFA)

HT-2 toxin (3,4-dihydroxy-15-acetoxy-8-[3-methyl-butyryloxy]-12,13-epoxytrichothec-9-ene (TFA).

Triacetoxyscirpene

8'-hydroxy-zearalenone (alpha isomer)

8'-hydroxy-zearalenone (beta isomer)

zearalenone (trans)

3'-hydroxy-zearalenone (F-5-2)

H-1 hemorrhagic factor produced by Fusarium (in the process of chemical characterization.

Fusarochromanone (TDP-1) underivatized (causes bone deformation)

Fusarochromanone (TMS ether)

TDP-2 TMS ether. Monoacetyl derivative of fusarochromanone

TDP-2 underivatized

TDP-3 unknown derivative of fusarochromanone(TDP-1)

**TDP-4 (Unidentified toxic derivative of Fusarochromanone)**  
**Monoacetyl derivative of Fusarochromanone(synthetic product)**

**B. USE OF FUSTOX LIBRARY FOR IDENTIFICATION(LIBRARY SEARCH-VG SYSTEM)**

An isolate of *Fusarium equiseti* which originated in Germany was cultured on rice with the intention of determining the presence of Fusarochromanone (Fusarium metabolite that causes tibial dyschondroplasia). A crude extract was obtained and resolved by capillary chromatography on the GC/MS VG-7070EQ. Figure 1 shows the resolution of the TMS derivatives on DB-5 capillary column. Multiple components were found (at least 5) that belonged to this chemical group. Scan 201 was suspected of being Fusarochromanone and was tested by the FUSTOX library. The results of the search of scan 201 is reported in figure 2. Note that among the hits or tentative identifications, entry number 8 is listed as the best fit which agrees with our own interpretation of the spectrum. The library is searched both forward (match the spectrum with the library) and reverse search (match the library with the spectrum). The higher the numbers in the table marked PUR, MIX and REV ( respectively 602, 833 and 654) the better the quality of the match. PUR stands for relative purity of the compound, i.e. very few flagged masses; MIX shows the results of the comparison of both the forward and reverse library search and REV shows the results of only the reverse search. A comparison of scan 201 with the library entry is shown in figure 3 ;the spectra are almost identical indicating an excellent implementation of the search system.

Scan #245 of figure 1 was also matched with the library and the results are presented in figure 4. Once more the library correctly identified the FUSAROCHROMANONE derivative (TDP-2). Note the high values of PUR, MIX and REV. The comparison of the library spectrum of TDP-2 with scan 245 is shown in figure 5; the spectra are almost identical. The search routine of the VG USER GENERATED LIBRARY appears to be excellent for implementation for *Fusarium* toxin identification.

The total ion current chromatogram of the extract from *Fusarium* is shown in figure 6. The TDP-1 indicator fragment is  $m/z+ 218$  and this ion was used to uncover naturally occurring derivatives of FUSAROCHROMANONE in the culture. At least 8 different derivatives were found. Of the eight, we only identified 2; the remainder are new components.






FIGURE 1. TOTAL ION CHROMATOGRAM (UPPER BOX) OF THE CRUDE EXTRACT FROM FUSARIUM EQUISETI OBTAINED FROM GERMANY. THE TMS ETHER DERIVATIVES WERE RESOLVED ON 10 METER CAPILLARY COLUMN (DB-5). (LOWER BOX): MASS SPECTRUM OF SCAN 201 INDICATING THE PRESENCE OF THE TMS ETHER OF FUSAROCHROMANONE (TDP-1).

LIB Library Search Report of DUTCHTDP#201 Library: FUSTOX 23-MAR-87

18 Library spectra compared for BEST FIT  
 18 matched 0 or more of the 16 largest peaks in the unknown

File	Entry	Compound Name	Reference
1	8	TDP-1 (TMS) FUSAROCHROMANONE-TMS 14OCT86 70ED 70EV EI	0-00-0
2	9	TDP-2(TMS DERIVATIVE) MONOACETYL DER OF TDP-1. 28AUG'86	0-00-0
3	7	TDP-1 (FUSAROCHROMANONE) 70ED EI 70EV TIC:4389900 DIP	0-00-0
4	13	MONOACETYL OF TDP-1 SYNTHESIZED BY XIEGEN. 28AUG'86 70E	0-00-0
5	10	MONOACETYL-TDP-1 OR TDP-2 UNDERIVATIZED 6AUG'86 70ED 7	0-00-0




FIGURE 2. REPORT OF A LIBRARY SEARCH MADE ON THE VG SYSTEM OF SCAN 201. NOTE THAT THE BEST FIT IS LISTED AS FUSAROCHROMANONE. THE PURITY IS REPORTED AS 602, RESULTS OF MIXED SEARCH (FORWARD AND REVERSE) IS 833 AND THE REVERSE SEARCH (REV) IS 654. THE HIGHER THE VALUE OF THESE NUMBERS THEN THE BETTER THE CORRELATION OR FIT. MAXIMUM VALUE IS 1000.

LIBFITS91: x1 Bqd=291 DUTCHTOP  
 TOP-1 (TRIS) FUSAROCROMANONE-TRIS 14OCT86 70EV 70E p682 R833 r654 RPN:8-00-0  
 LIB: FUSTOX 8 Bpk: 347 Mt: 347




FIGURE 3. COMPARISON OF THE MASS SPECTRA OF THE SCAN 201 (UPPER BOX) WITH THE LIBRARY ENTRY OF FUSAROCROMANONE. THE SPECTRA ARE ALMOST IDENTICAL INDICATING AN EXCELLENT FIT.

LIB Library Search Report of DUTCHTOP245 Library: FUSTOX 23-APR-87

18 Library spectra compared for BEST FIT  
 18 matched 8 or more of the 16 largest peaks in the unknown

File	Entry	Compound Name	Reference
1	9	TDP-2(TRIS DERIVATIVE) MONOACETYL DER OF TDP-1. 28AUG'86	8-00-0
2	8	TOP-1 (TRIS) FUSAROCROMANONE-TRIS 14OCT86 70EV 70E	8-00-0
3	14	TOP-4 (DIACETYL-TDP-1) SYNTHESIZED BY XIEGH. 3 JUN'86	8-00-0
4	7	TOP-1 (FUSAROCROMANONE) 70EV EI 70EV TIC84393000 DIP	8-00-0
5	13	MONOACETYL OF TDP-1 SYNTHESIZED BY XIEGH. 28AUG'86 70E	8-00-0



FIGURE 4. MASS SPECTRAL LIBRARY SEARCH OF THE E.I. SPECTRUM OF SCAN 6245 TENTATIVELY SUSPECTED OF BEING TDP-2 (MONOACETYL DERIVATIVE OF FUSAROCROMANONE). TDP-2 WAS SELECTED CORRECTLY AS THE BEST FIT ; THE PURITY, MIXED SEARCH AND REVERSE SEARCH VALUES WERE HIGH INDICATING A GOOD MATCH.

LIBRARY#1: x1 800-245 DUTCHTOP  
 TOP-2CTRS DERIVATIVE MONOACETYL DER OF TOP-1. 20 p614 10790 r725 RFN:0-00-0  
 LUB:FUSTOX 9 Bpk: 218 MA: 406




FIGURE 5. COMPARISON OF THE MASS SPECTRA OF SCAN 8245 OF THE CRUDE EXTRACT AND THE LIBRARY COPY OF THE TOP-2 STANDARD. CASUAL INSPECTION WOULD INDICATE THAT THE SPECTRA ARE ALMOST IDENTICAL.

DUTCHTOP 01-321 23-APR-87 17:25 7000 Port:  
 A: RTIC B1: 218 C1: 252 D1: 354 E1: 347 F1: 406 G1: 244 H1: 411  
 Tack: GERMAN ISOLATE FOR TOP-1 AND TOP-2 1802R TFS DERV. GC/00/30

Sub: PARENT




FIGURE 6. THE TOTAL ION CHROMATOGRAM OF THE CRUDE EXTRACT OF THE FUSARIUM CULTURE IS SHOWN IN BOX A. BOX B SHOWS THE TOTAL ION CHROMATOGRAM OF MASS 218 WHICH IS THE BASE PEAK OF TOP-1 AND AN EXCELLENT INDICATOR OF THE TOP-1 NUCLEUS. THE TRACING IN BOX A INDICATES THAT AT LEAST 7 DIFFERENT DERIVATIVES OF FUSARICHRONAMONE ARE PRESENT IN THIS MIXTURE. SCANS 186, 201, 202, AND 245 ARE MOST ABUNDANT. SCANS 201 AND 245 ARE TOP-1 AND TOP-2 RESPECTIVELY; THE OTHER COMPONENTS HAVE NOT BEEN IDENTIFIED.

## INTERFERING SUBSTANCE IN HUMAN BLOOD AND ITS ANALYSIS BY THE DAUGHTER ION MASS SPECTRAL LIBRARY

During the course of analysis of human urine for various metabolites of T-2 toxin, it came to our attention that a substance present in control blood was interfering with our analysis of T-2 toxin by Multiple Reaction Monitoring (MRM). The material did not actually interfere with T-2 because it had a different retention time on the capillary column but theoretically it was not supposed to pass through the magnet. By way of background, in MRM analysis of T-2(TFA), the magnet is parked on  $m/z+478$  and the quadrupole is scanning for selected daughter ions (the quadrupole is programmed to detect only the legitimate daughters of the compound in question). This means that either the daughters are there or they are not; the quadrupole will not detect alien daughter fragments. The magnet end is set on the accurate mass of 478.14506 to insure that it sits directly on top of the selected parent; this does not mean that because of accurate mass selection that the instrument is scanning at high resolution; it is not. Normally we operate at about a resolution of 1000. One does not expect any appreciable noise to filter through the system. However, in the analysis of blood, another 478 mass fragment appeared at a retention time of about 30 seconds past that of T-2 (figure 7-A). The peak was appreciable and at low concentrations of T-2, caused the T-2 peak to look insignificant because of its normalization to the most abundant peak. We decided to analyze the interfering peak with the daughter ion library and to compare it with daughters of T-2-TFA.

The total ion chromatogram of the interfering substance is shown in figure 7. The peak is shown in scan 134 and its daughter ions are shown in the lower box. We could not detect any  $m/z+478$  fragment in the daughter ion scan and could not detect any fragments characteristic of T-2 -TFA daughters. As shown in figure 8, the DAUGHTER ion mass spectral library did not generate any match for the interfering blood substance. In contrast, a sample containing T-2 toxin was run in the daughters system and its daughter ion spectrum was matched with the library. As noted in figure 9, a good fit was found in the library entry and all the daughter ions were found. A comparison of the library entry with that of the unknown is shown in figure 10. Note the identical daughter ion spectrum of mass 478 of the library and the sample. The parameters of the library scan i.e. values used in directing the search, are shown in figure 11. It is important to adjust the fit threshold (THR) to 100 or lower. It is also helpful to enable the filters to select only those samples that have the correct daughters.




FIGURE 7. DAUGHTER ION SPECTRUM OF AN INTERFERING SUBSTANCE PRESENT IN CONTROL HUMAN BLOOD WITH A  $M/E$  478 FRAGMENT SIMILAR TO THAT OF T-2-TFA. NOTE THAT THE 478 FRAGMENT IS RAPIDLY DECOMPOSED AS NO TRACE OF IT IS FOUND ONCE IT LEFT THE COLLISION CHAMBER (CAD) AND THAT THE DAUGHTERS ARE DIFFERENT FROM THOSE FOUND IN T-2-TFA. THE MAGNET WAS PARKED ON MASS 478.145 AND THE QUADRUPOLE WAS SCANNING BETWEEN MASSES 700 AND 50. THE RETENTION TIME OF T-2 TFA IS ABOUT 6.2 (EQUIVALENT TO SCAN 117) AND ITS ANALYSIS IS NOT INFLUENCED BY THE UNKNOWN SUBSTANCE SHOWN ABOVE. HOWEVER, THE SUBSTANCE IS DESCRIBED BECAUSE CONTRARY TO EXPECTATION, IT HAS AN IDENTICAL MASS OF 478.

2PPB312 12-MAR-87 16:57 70EQ Acnt: SUS:MKCM  
 GR 1 A: 478.1450 B: 478.1450 C: 478.1450 D: 478.1450  
 Text: SPIKED BLOOD 2PPB T-2 TFA DERV. 478-478,180,138,121




FIGURE 7A. Resolution of the T-2 toxin TPA from the interfering substance found at retention time 7.0 whereas T-2-TFA is found at 6.25. Note that the T-2 concentration is equivalent to 200 picograms which in box A (fragment 478 is monitored) is actually in a concentration higher than T-2. However, boxes B, C and D are true daughters of T-2 (180, 138 and 121 respectively) and no problem was encountered in detection.

LIB Library Search Report of BLDSUB#134

Library: DAUGHTER

13-APR-87

6 Library spectra compared for BEST FIT

6 matched 0 or more of the 16 largest peaks in the unknown

Flt Entry Compound Name  
 NO FITS FOUND

Reference

Flt Elements Bpk Mt pur MIX rev




FIGURE 8. RESULTS OF A LIBRARY SEARCH OF THE DAUGHTER ION MASS SPECTRAL LIBRARY OF TRICHOETHENES. NO FIT OR MATCH WAS FOUND BECAUSE THE DAUGHTERS OF THE UNKNOWN DID NOT CORRESPOND TO THOSE OF T-2-TFA. THIS SEARCH SERVES AS A TEST TO PROVE THAT ALTHOUGH AN UNKNOWN MAY HAVE A FRAGMENT SIMILAR TO T-2-TFA, IT CAN BE EASILY BE DISCRIMINATED FROM T-2-TFA BY MEANS OF CHARACTERISTIC FRAGMENT IONS. F

LIB Library Search Report of 478T2876

Library: DAUGHTER

13-APR-87

6 Library spectra compared for BEST FIT

6 matched 0 or more of the 16 largest peaks in the unknown

Flt	Entry	Compound Name	Reference
1	6	DAUGHTERS OF 478 OF TFR T-2. 70EQ 70EV EI TIC 684a 200	0-00-0
2	1	DAUGHTERS OF 400 OF T2-TFR. 29JUL '86 70EQ 70EV EI TICS	0-00-0

Flt	Elements	Bpk	Mt	pur	MIX	rev
1		121	0	360	496	546
2		400	0	12	97	101




FIGURE 9. RESULTS OF THE LIBRARY SEARCH OF THE DAUGHTER ION MASS SPECTRAL LIBRARY FOR T-2-TFA. THE LATTER WAS DETECTED IN A URINE SAMPLE AND ITS IDENTITY VERIFIED BY MEANS OF THE LIBRARY. NOTE THE HIGH VALUES IN THE MIX AND REV PARAMETERS WHICH INDICATE A HIGH PERCENTAGE PROBABILITY OF IDENTIFICATION OF THE T-2-TFA.

LIBFITS#1\* x1 Bgd=76 478T2 p388 M496 r546 RFN:0-00-0  
 DAUGHTERS OF 478 OF TFA T-2. 70EQ 70EV EI TIC 682 LLb DAUGHTER 5 Bpk 121 Mut. 0




FIGURE 10. COMPARISON OF THE DAUGHTER ION LIBRARY SPECTRUM OF T-2-TFA WITH THE UNKNOWN FOUND IN THE URINE SAMPLE SHOWN IN FIGURE 9. THE COMMAND "SN" ACTIVATES THIS COMPARISON. NOTE THE IDENTICAL MASSES IN EACH. THE INTENSITIES OF FRAGMENTS 478, 180 AND 121 VARY BECAUSE OF THE DIFFERENCE IN COLLISION ENERGY USED IN EACH DETERMINATION. IT IS ACCEPTED THAT COLLISION ENERGY WILL VARY AND THUS THE ABUNDANCE OF DAUGHTERS IONS WILL BE AFFECTED; HOWEVER, THIS DOES NOT APPEAR TO BE SERIOUS ENOUGH TO INFLUENCE THE LIBRARY SEARCH.

LIB Library Parameters  
 MS Mass range 40, 1000  
 MS Peaks to consider 16  
 MS Machine data Y  
 MS Peaks per 100 ms 40  
 MS Pass 1 window 50  
 MS Pass 2 window 7  
 FIT Search type (P,B,D) A  
 THR Filz threshold 50  
 DF Min. in data search 20  
 DF Excl. of masses  
 FL Full library Y  
 LL Library limits 6, 40000  
 BGD Background edit mode A  
 LIB List weighted ions N  
 UPS Keep pre-ach /S  
 PR Manual reporting Y  
 PR Manual report /R  
 UPS Keep pre-ach /R  
 PR Manual p.c.m. /R

LIB Library Parameters Page 2 of 2  
 GENERAL FILTERS  
 OF Mass filters E  
 MC Molecular weight 0, 0  
 E1 C12,H14,N2,O5 range C:0.511 N:0.122 O:0.53 S:0.53 Cl:0.53  
 E2 C12,H14,N2,O5 range F:0.53 P:0.53 S:0.53 Cl:0.53  
 E3 C12,H14,N2,O5 range  
 DF Mass filter  
 DF Isomer filter 0  
 CF Classification flag 0  
 EDIT FILTERS  
 MS Compound name  
 MS Excluded masses 121, 180, 138  
 MS Reference number

FIGURE 11. PARAMETERS USED IN THE LIBRARY SEARCH ROUTINES WHICH ALLOW VERSATILITY IN A SEARCH. AS AN EXAMPLE, THE PARAMETER "ROM" CAN BE ADJUSTED SO THAT REQUIRED DAUGHTERS MUST BE PRESENT IN THE SPECTRUM BEFORE THE SEARCH ROUTINE WILL CONSIDER A MATCH. IN THIS EXAMPLE, M/Z+ 121, 180, AND 138 WERE SELECTED AS MASS FILTERS. THE PARAMETER (DF) MUST BE ENABLED (ENABLE FILZPT) BEFORE DISCRIMINATION WILL BE ACTIVATED. ACCORDINGLY THE FILTERS CAN BE DISABLED BY TYPING IN "D". THE FIT THRESHOLD "THR" TIGHTENS THE SEARCH PARAMETERS BY CAUSING A GREATER DEGREE OF SELECTIVITY AS ITS VALUE IS DECREASED. IN THIS EXAMPLE, A THR VALUE OF 50 WAS SELECTED.



## ANALYSIS OF T-2-TFA USING MULTIPLE REACTION MONITORING

The use of multiple reaction monitoring (MRM) for trace analysis of T-2 toxin is based on prior experimentation in which appropriate parent ions are studied and selected for use in the analytical system. As an example, the daughters of T-2-TFA were determined and then selected for ultimate sensitivity in the assay. Thus, as one examines the full spectrum of T-2-TFA in E.I. (figure 12) it suggests that certain fragments be studied. The molecular ion (562) is normally too weak to be seen but after magnification (5X), it is readily visible. Fragments such as  $m/z+461$ , 478, 401 and 327 appear significant enough to yield acceptable analytical results. Figure 13 shows the daughter ions of  $m/z+401$  and careful inspection indicates that although the parent ion is appreciable, the daughters are very weak. In contrast, the daughters of  $m/z+478$  give very intense fragments and have all the characteristics of a good ion to monitor i.e. strong daughters in the low mass region and low signal to noise ratio (figure 14). The fragments of interest are: 478, 180, 138 and 121. These are then used qualitatively and quantitatively in subsequent tests for the analysis of trace (2ppb) of T-2-TFA in the biological matrix of blood.

Human blood samples were amended with the equivalent of 2ppb of T-2 toxin and extracted and cleaned by our general laboratory procedure. The extract contained a biological matrix; the trifluoroacetate derivative of T-2 was made in situ and injected into the VG GC/MS operating in the MRM mode. The calculated amount injected into the capillary column was approximately 160 picograms after accounting for a 20% loss in extraction. Note the excellent sensitivity and base line when the  $m/z+180$  daughter is monitored (figure 15). The analysis and quantitation of daughters 180, 478, 138 and 121 are shown in figure 16. The multiple reaction monitoring system of analysis appears to be the best yet developed in our laboratory for T-2 toxin analysis. The most important determinant in the analysis is the lack of spurious peaks and electronic noise when working at this level. Our signal to noise ratio at 2 parts per billion is about 100:1, the best we have ever achieved. We recommend this method of analysis because of its sensitivity and unequivocal identification of T-2 toxin.

STDS#216 x1 Bgd=212 28-MAR-87 19:31-0:05:17 70EQ EI-  
 BpM=0 I=1.2v H=0 TIC=82454000 Acnt: Sys: PARENT  
 Text: STDS T2, HT2, T240L, TC1, TC3, CTFD 50NG/UL, 28MAR '87 PT=0 Cal: ACCU2  
 #216 5.0  
 8094000  
 1998000




FIGURE 12. FULL SCAN MASS SPECTRUM OF T-2 -TFA IN E.I. AT 70 EV AND A SOURCE TEMPERATURE OF 150 C. NOTE THE MOLECULAR ION (562) WHICH NORMALLY IS NOT VISIBLE. WE USE A SOURCE TEMPERATURE OF 200 TO ACHIEVE THESE REPRODUCABLE RESULTS. THE FRAGMENTS OF INTEREST ARE 478, 460, 401 AND 327.






FIGURE 13. DAUGHTER IONS OF FRAGMENT 401 . NOTICE THAT VERY FEW FRAGMENTS OF INTEREST ARE GENERATED AND ALL ARE OF LOW INTENSITY.

2PPB312 12-MAR-87 16:57 70EQ Acnt: Sus:ARM  
 GR 1 A: 478.1450 B: 478.1450 C: 478.1450 D: 478.1450  
 Text: SPIKED BLOOD 2PPB T-2 TFR DERV. 478-478.180.138.121



T2478A#1 x1 F:D FFR:1 29-JUL-86 01:28-0:00:00 70EQ EI+  
 BpA=0 I=57mv Ha=700 TIC=1495000 RV Acnt: Sys:DAUGHTERR  
 Text: Average of DR??:T2478#18-24 PT=0 FR:0.000000



2PPB312 12-MAR-87 16:57 70EQ Acnt: SUS:NRN  
 GR 1 A: 478.1450 B: 478.1450 C: 478.1450 D: 478.1450  
 Text: SPIKED BLOOD 2PPB T-2 TFA DERV. 478-478,180,138,121




FIGURE 16. ANALYSIS OF 2PPB OF T-2 TFA IN HUMAN BLOOD USING MULTIPLE REACTION MONITORING. THE PARENT FRAGMENT USED IS 478 AND THE DAUGHTERS MONITORED ARE 478, 180, 138 AND 121.

#### DAUGHTER IONS OF THE TFA DERIVATIVES OF T-2, HT-2 AND T-2-TETRAOL


Monkey urine samples were received during the first quarter of 1987 with the intention of analyzing for T-2 toxin and its derivatives during various times after dosage of the animals. In order to accomplish this task, it was necessary to first determine the daughter ions of some of the expected derivatives in order to analyze by multiple reaction monitoring. This task has been accomplished for T-2, HT-2 and T2-tetraol. The daughters are shown in the following figures.

Figure 17 shows the full mass spectrum of T-2-TFA and the daughters of the major ions. Note that the molecular ion 562 gives rise to 478, 460 and 461. The daughter 400 is reached by multiple reactions from 460 and 401. The parent of choice in the MRM analysis is 478.


Figure 18 shows the full mass spectrum in E.I. of HT-2 toxin (TFA). The molecular ion 616 fragments to 532 which in turns gives rise to 180, 138, and 121 as in T-2 except that the origin of these daughters in T-2 is from 478. The parent ion 532 is the fragment of choice for MRM analysis.

Figure 19 shows the full mass spectrum of the TFA derivative of T-2 -tetraol. The molecular ion is seen in E.I. but is not very intense although if detected it gives rise to very intense daughters. Parents 568 and 569 give rise respectively to daughters 454 and 455. Parent 330 is very intense and gives rise an intense 217 daughter and is the choice in the analysis of T-2-tetraol by MRM.


STD38216 x1 Bgd=212 28-APR-87 19:31-0:05:17 78C3 EI-  
 Sp=0 I=1.2v Hs=0 TIC=82454000 Rank: Sys: PARENT  
 Text: STDS T2, HT2, T24OL, TC1, TC3, CTFD58NG/UL, 28APR 87 PT=0 Cal: ACCUR  
 8216 1.0  
 8894800  
 1393000



STD38183 x1 Bgd=177 28-APR-87 19:31-0:04:31 78C3 EI-  
 Sp=0 I=2.2v Hs=0 TIC=83817000 Rank: Sys: PARENT  
 Text: STDS T2, HT2, T24OL, TC1, TC3, CTFD58NG/UL, 28APR 87 PT=0 Cal: ACCUR  
 8182 1.0  
 1788000



STD38188 x1 Bgd=186 28-APR-87 19:31-0:02:47 78C3 EI-  
 Sp=0 I=2.3v Hs=0 TIC=146817000 Rank: Sys: PARENT  
 Text: STDS T2, HT2, T24OL, TC1, TC3, CTFD58NG/UL, 28APR 87 PT=0 Cal: ACCUR  
 8188 1.0  
 1538000



# ANALYSIS OF T-2 TOXIN IN HUMAN BLOOD USING TANDEM MASS SPECTROMETRY GAS CHROMATOGRAPHY

## INTRODUCTION

The utilization of combined gas chromatography tandem mass spectrometry for the analysis of T-2 toxin in blood is presented here. The analysis was carried out to confirm and quantitate the presence of T-2 toxin at a low part per billion level (2ppb in 2mL blood). The procedure involved spiking fresh whole blood with the toxin, extraction with an organic solvent, derivatization with trifluoroacetic acid anhydride, and GC/MS/MS analysis. Please remember that the following work was carried out in a biological matrix consisting of the equivalent of 2mL of human blood.

## AMENDING BLOOD WITH T-2 AND EXTRACTION

Blood samples (2mL) were amended with T-2 (1ng/ul) in methanol. They were centrifuged for 2 minutes at 7000 rpm in 30 mL tubes with 0.5 mL KH<sub>2</sub>PO<sub>4</sub> and 10 mL of toluene:acetonitrile (95:5 v:v). The organic layer was pipetted into a round bottom flask and evaporated on a rotary evaporator. Samples were transferred to half dram vials with methanol and evaporated under nitrogen with gentle heat.

## DERIVATIZATION

The dried down samples were derivatized with trifluoroacetic acid anhydride (100 uL) and heated to 60 degrees C. for 20 minutes. After which the excess reagent was evaporated under nitrogen. Samples were brought up in 20 uL of toluene for injection. 1 ul was taken for injection.

## GC CONDITIONS

Samples were injected onto a 10 meter fused silica capillary column via a Hewlett Packard on column injection device. The injection device is shown in figure 1. The oven temperature was programmed from 80 to 300 at 30 deg. per second. T-2 eluted at 6.28 minutes into the run. The helium gas flow was set to 10 psi.

In addition to the on column injection technique a 2 ppb blood sample was prepared for splitless injection. The technique afforded overall better resolution and narrowly defined peaks. Signal to noise ratios for 2ppb samples are better than 14 to 1. This method proved superior to the present on column technique. Injection was done on a 10 meter capillary column. Column purge was activated at 0.6 minutes after injection.

## SOURCE CONDITIONS

A standard VG 7070 EI/CI source was used in the analysis. The instrument was operated at 6kV, 70 electron volts. Source temperature was set at 100 degrees. Resolution was set to 2000.

#### MS/MS TUNING

The instrument was tuned for MS/MS analysis with perfluorotributyl amine. Intermediate and quad lens settings were set to allow maximum transmission of parent ion 614. The quadrupole mass analyzer was calibrated with the VG software LOW MASS (LM) and HIGH MASS (HM) ion commands. The high resolution portion of the mass spectrometer was tuned with perfluorokerosene using the CAL programs.

#### COLLISION ACTIVATED DECOMPOSITION (CAD) set-up:

Argon was used as the collision gas and the chamber was pressurized to  $1 \times 10^{-6}$  torr. Collision energy was set at 15 volts.

#### FA3 AND MULTIPLIER SET :

The gain was set on FA3 unit to  $2 \times 10^{-7}$ . Response time was set to 0.01 and multiplier set to 2.10.

#### MULTIPLE REACTION MONITORING PROGRAM (MRM):

Selected Ion Recording (SIR) was called up and MRM was set up to monitor the decomposition products (daughters) of ion 478 from T-2 TFA, to 180, 138 and 121 as seen in the figures. The exact mass 478.145 was used to set the parent ion correctly. When running the program we found it better to set the magnet manually by positioning it over the parent ion rather than letting the data system do this. To do this one must allow the data system to have control over quad scanning only. The 7070EQ is equipped with knobs on the scan control unit to select this type of scanning procedure. The parameter settings are outlined in figure 4.

#### RESULTS:

The electron impact mass spectrum is shown in figure 2. Please note that the molecular ion ( $M^+=562$ ) is visible in the spectrum. This was accomplished by lowering the source temperature to 100 degrees C. This is not the normal procedure in EI as the normal operating temperature is 200 C. Our source heater was not operative in this situation and this allowed us to reach the lower temperature setting. We are not recommending this temperature setting. A typical ON-COLUMN injection resolution is found in figure 5. The ions monitored were 180, 138 and 121 and their ratios as daughters are shown. The ratios were satisfactory however the resolution of the components of the standard did not meet out criteria of acceptability. The peak sensitivity and peak width was not acceptable. Figure 5 shows resolution of T-2 TFA as the pure



standard whereas figure 6 shows T-2-TFA in a biological matrix.

The parent ion 478 gave distinct and intense daughter ions that (figure 3) allowed elimination of other components from the background except for one present in human blood and having the exact mass (478) as T-2-TFA toxin. This metabolite (figures 7 and 8) found in these blood samples did not interfere with the analysis of T-2 because it resolved itself nearly one minute from T-2. In addition, the ion ratios for the interference did not fit the T-2 toxin and therefore would not be confused with T-2 (figure 7). However, it is mysterious that this fragment should pass through the magnet as the resolving power was set to three decimal places. The major problem in the analysis of 1 or 2 ppb quantities at this time is quantification and reproducibility. One part per billion is detectable based on the 2 part per billion spiked sample which had a 20 to 1 signal to noise ratio (figure 8). The ratio variance for samples and standard injections was only 10%. This means that the ratio of ions can be used to confirm and quantitate the toxin. With such a small sample size (2 gram blood) and low amounts injected (400 picograms) reproducibility suffered much. Problems may be due as well to sample injection (use of the syringe as well since the matrix was rather viscous). An internal deuterated standard or homolog would be necessary for quantifying these amounts of toxin.

On column chromatography had some benefits and draw backs in the analysis. It allowed for overall greatest sample throughput into the mass spectrometer but peak broadening, column coating degradation, and increased tailing may severely limit its performance on particularly dirty samples. The SPLITLESS method of sample introduction (figures 8 through 12) gave excellent chromatography for analysis in this low concentration region. Interpretation of data: The data appeared to be as easy to interpret as that of normal selected ion analysis. The ratio consistency as mentioned above was quite satisfactory. The purpose in using this method of analysis (MRM) was to demonstrate the lack of interferences and sample limits of detection. Base lines such as those shown in these analyses are quite uncommon for mixed matrix analyses. We were quite pleased to eliminate other chemical noise from the base line. This procedure should also be compared to other such procedures as selected ion recording.


The splitless method of injection is demonstrated in figure 9 which shows the resolution of 2ppb of T-2-TFA in a biological matrix of the equivalent of 2ml of blood. The daughter monitored was  $m/z + 180$  which showed a minimum of noise ( $S/N = 17:1$ ). The daughters 138 and 121 are shown in figures 10 and 12 respectively. Their signal to noise ratios (12:1 and 14:1) are acceptable but inferior to that of 180. Our recommendation is that all three daughters be monitored for the purposes of calculating ratios but that 180 be used for quantitation.

Figure 12 shows the resolution of T-2-TFA in MRM using the

splitless method of injection. Injection of one nanogram under ideal conditions shows excellent resolution and peak width superior to that found in the On-Column injection shown in figure 5.

We conclude from these studies that the SPLITLESS method of sample introduction into the GC/MS is the method of choice. It gives superior resolution when compared to ON-COLUMN injection. Mass fragment 478 should be used as the parent ion for analysis and the daughter 180 gives the best quantitative results. Fragment 478 was chosen after careful study of the daughters of 460, 461, 401 and 418. All of the above was done in the electron impact mode of analysis. The latter is traditionally less sensitive than positive chemical ionization (PCI) in methane which is our next effort in this study. We expect to increase our sensitivity by a factor of ten due to the large fragment ions encountered in PCI.

FIGURE 1. On-Column injection device used in the injection of the T-2 amended blood matrix through the Hewlett Packard gas chromatograph. This device consists of a 0.53 mm I.D. precolumn (about 15 cm in length) coupled to 0.25 mm I.D. 10 meter DB-5 analytical column. This method allows direct injection with a regular metal injection syringe rather than a flexible on-column injection syringe. The device is patterned on that of Hartman et al. of Rutgers University who has developed the procedure followed.



STAN86 x1 Bgd=22 10-FEB-87 15:30-0:07:12 70E0  
 I=423eV H=0 TIC=22308000 Acnt:  
 st:T-2 TFA ON-COLUMN

EI-  
 Sys: PARENT  
 PT=0 Cal: ACCU1

#86 1.0  
 2816000  
 55000





FIGURE 3. Daughters ions obtained from fragmentation of ion 478 from the TFA derivative of T-2 toxin. Ions 180, 138 and 121 were chosen for subsequent MRN analysis. The fragment 478 was found to be superior (gave less background noise) to that of 401 in multiple reaction monitoring as determined in previous analyses.



* SYSTEM:MRM		Parameters for Selective Ion Recording [G# 1 MRM] (Sector)			
DAT	Data filename	C:TESTB	IAV	Maximum volts	6000
REF	Reference filename	PFK	IMR	Maximum mass at IAV	2256
INS	Instrument	1:70EQ	GTM	Time 0:06:00 0:20:00	Mode EI+B Gas
ACH	Customer account		CHN	M(amu) S(ms) D(ms)	M(amu) S(ms) D(ms)
ACV	Accelerating volts	6000	478	180	100 10
RES	Instrument resolution	1000	478	138	100 10
SAP	#Samples	1	478	121	100 10
INJ	#Injections	1			
GRP	#Groups	1			
CLS	#Calibration scans	0			
CST	Cal. scan time (s)	5			
CTL	Cal. tolerance (ppm)	200			
CEX	Cal. examination	Y			
PEX	Peak examination	Y			
LMS	Lock span (peak widths)	2.0			
LST	Lock step (peak widths)	0.02			
FLO	Fast Lock on	N			
FIGURE 4. Multiple Reaction Monitoring Parameter Page For T-2 toxin TFA Derivative. The exact parent ion mass of 478.145 was selected. The daughters monitored were 180,138, and 121. .PA					
TXT	Sample# 1: 500 PICO T-2 TFA DEV. MRM ON-COLUMN				
'H=hardcopy RETURN=next ESC=prev CTRL/A=abort					
'G=go 'Q=quit 'C=create 'DEL=delete 'O=overwrite 'Z=zero <group,sample>					

ESTB 11-FEB-87 21:40 70EQ Acnt:  
 R 1 A: 478.1450 B: 478.1450 C: 478.1450  
 ext: T-2 TFA MRM ON-COLUMN

Sus:MRM




FIGURE 5. Ratios and retention time of T-2 TFA pure standard. Retention time 6.28 minutes. Ions 121/180 0.46, ions 138/180 0.53 taken from parent ion 478 going to ions 121, 138, and 180. On-column injection. 10 meter 0.25 capillary column. GC oven temperature program 60/100 at 10 deg. per sec.

TESTB 12-FEB-87 00:13 70EU Hunt: 545:MKM  
 GR 1 A: 478.1450 B: 478.1450 C: 478.1450  
 Text T-2 SPIKED BLOOD 2N/G 2G EQUIVALENT TFA DERV. IN 20UL




FIGURE 6. Multiple reaction monitoring (MRM) of a 2ul blood sample amended with 2ppb of T-2 toxin. The sample was introduced into the GC/MS via ON-COLUMN INJECTION. This method of injection will be contrasted with splitless injection in a later section.

GR 1 H: 478.1450 B: 478.1450 C: 478.1450  
 Text: ~~XX~~




FIGURE 7. Component found in donor blood having similar properties to that of T-2 toxin except that its retention time is different (Rt=4.15 for T-2 vs 4.1 for the unknown). The daughters of the unknown have different ratios from that of T-2 (T-2= 121/180=0.45 and 138/180=0.75; ratios of unknown= 121/180=1.34 and 138/180=0.74). The unknown component has a significant m/z=478 similar to T-2-TFA so that when the magnet is parked on 478, both the T-2 and the unknown component are passed through.

2PPBT23 17-FEB-87 17:45 70EQ Acnt: Sus:NRN  
 GR 1 A: 478.1450 B: 478.1450 C: 478.1450 D: 478.1450  
 Text: S/N 14.2/1  
 100 C 1, 1, 1




FIGURE 8. Resolution and quantitation (2ppb in blood matrix) of T-2 TFA toxin by MRM using the SPLITLESS method of compound introduction. Note the 10 second peak width at the base which is a great improvement over the ON-COLUMN injection. The unknown metabolite in the blood sample can be seen at Rt 8.3; its resolution has also been increased. Compare with resolution as shown in figures 6 and 7.

2PPBT23 17-FEB-87 17:45 70EQ Acnt: Sus:NRN  
 GR 1 A: 478.1450 B: 478.1450 C: 478.1450 D: 478.1450  
 Text: S/N 14.2/1






FIGURE 9. Resolution and quantitation of 2ppb T-2 toxin by MRM using the SPLITLESS method of component introduction. The daughter ion monitored was m/z 180; its signal to noise ratio is about 17:1. Fragment 180 appears to be the most sensitive for monitoring T-2-TFA.

2PPBT23 17-FEB-87 17:45 70EQ Acnt: Sus:NRN  
 GR 1 A: 478.1450 B: 478.1450 C: 478.1450 D: 478.1450  
 Text: S/N 14.2/1



2PPBT23 17-FEB-87 17:45 70EQ Acnt: Sus:NRN  
 GR 1 A: 478.1450 B: 478.1450 C: 478.1450 D: 478.1450  
 Text: S/N 14.2/1





TESTB 17-FEB-87 15:26 78EB Acnt: Sys:NRN  
GR 1 A: 478.1450 B: 478.1450 C: 478.1450 D: 478.1450  
Text: T-2 TFA DERV. NRN EXP 478-180.138.121 1NG/UL  
100% A 1, 1, 1




FIGURE 12. Standard T-2 (1ng/ul), splitless injection using NRM as the analytical mode.

**TOXICITY AND TOXINS PRODUCED BY ISOLATES OF FUSARIUM  
SPECIES OBTAINED FROM NORWAY AND NEW ZEALAND**

Table 1 shows the complete record of the toxicity of the isolates of *Fusarium* obtained from Norway in 1984. The collections were made from the Arctic as well as the central and southern regions of that country. The species most prevalent were *Fusarium acuminatum*, *F. avenaceum*, *F. oxysporum* and *F. sambucinum*. From these groups, *F. oxysporum* and *F. sambucinum* emerged as the most toxic although toxicity was also found among fewer isolates of *F. acuminatum* and *F. avenaceum*.

All of the toxic isolates shown in table 1 were saved and examined for the kinds and quantities of known toxins they produce in order to account for the toxic lesions. The toxins that were analyzed for are: trichothecenes; HM-8, a cytotoxic unknown *Fusarium* metabolite presently being characterized in our laboratory; moniliformin; Fusarin-C; and H-1 the toxic hemorrhagic agent presently being characterized in our laboratory.

Table 2 shows the various known toxins produced by the various isolates of *Fusarium* in an effort to account for their toxicity. Within the *F. acuminatum* group, the predominant toxins found were HM-8 and moniliformin. *F. avenaceum* produced HM-8, moniliformin and infrequently Fusarin-C; only one isolate produced the hemorrhagic agent known as H-1. Although a considerable portion of the isolates of *F. oxysporum* yielded HM-8 and moniliformin, the majority of the toxicity could be accounted for by H-1. It is curious to note that none of the isolates that produced H-1 produced HM-8 as well and vice versa. Within the taxonomic group known as *F. sambucinum*, the largest percentage produced H-1. Only one isolate produced HM-8 and H-1 together.

Table 3 shows a compilation of the various toxins that were found in the toxic groups of *Fusarium*. In these analyses, TDP-1 (fusarochromanone) was also looked for. TDP-1 is the *Fusarium* metabolite that causes bone disorders. None of the isolates investigated produced TDP-1 whereas the majority of the *F. sambucinum* isolates produced H-1. The *F. oxysporum* group competed well with *F. sambucinum* for the production of H-1; however, the members of this group also produced moniliformin but not simultaneously with H-1. The significance of this observation is not apparent at present. The production of moniliformin was well represented among *F. acuminatum* and *F. avenaceum*.

An effort was made to determine whether the H-1 producing isolates could be found in regions outside of Norway. A collection was made in Tara Hills, a mountainous region in the South Island of New Zealand. Six isolates, all of them belonging

to the *Fusarium sambucinum* group, produced this toxin (Table 4). The Tara Hills are characterized as high grassy hills some of which receive a lot of shade; the climate is cold in the winter. The sheep grazing in this area become unthrifty and often called "poor doers" and this was particularly associated with the shady region of the Tara Hills complex. It is possible that *Fusarium* colonizing the grass may influence the health of these animals by the production of toxins such as H-1.

The isolates examined in Table 4 were selected because of their toxicity to rats, our bioassay test animal. Toxicity was accounted for by those isolates (6) that produced H-1, four isolates that produced moniliformin, and the six isolates that produced fusarenone-x. The toxins found ranged from deoxynivalenol to 3-acetyldeoxynivalenol, HM-8, zearalenone, moniliformin and H-1. Isolate 6F4 is being studied to identify the toxic component accounting for death in the test animals. This isolate produced DON, 3-acetyl-DON and zearalenone; these do not account for the toxicity.

Table 1. Toxicity of Fusarium isolates in rats fed a 1:1 mixture of Fusarium culture (rice) and complete rat diet.

<u>Fusarium</u> species	# isolates	# isolates causing					Food consumption (g) <sup>c/</sup>
		Death	Hemorrhage <sup>a/</sup>	Wt. loss (-)	No change (0)	Wt. gain (+)	
Control 1 (n=30)	-	0	0	0	0	30	43.13 ± 3.01
Control 2 (n=30)	-	0	0	0	0	21	41.33 ± 1.10
<u>F. acuminatum</u>	25	4	0	23	0	2	23.22 ± 3.13
<u>F. avenaceum</u>	28	7	3	20	1	7	20.51 ± 3.08
<u>F. culmorum</u> <sup>b/</sup>	1	0	0	1	0	0	17.30 ± 3.64
<u>F. oxysporum</u>	31	25	19	28	0	3	10.55 ± 2.74
<u>F. sambucinum</u>	40	37	33	40	0	0	7.88 ± 1.46
Total # isolates	125	73	55	112	1	12	

<sup>a/</sup> Includes hemorrhage of stomach, intestine, thymus, as well as hematuria.

<sup>b/</sup> This isolate caused uterine enlargement.

<sup>c/</sup> Each value is the mean for the total number of rats used for each species ± standard deviation.

Table 2. Production of mycelium from *Fusarium* isolates grown on rice medium.

<i>Fusarium</i> species and code number	Source	Area of origin	Mycelium*
Control 1	Complete rice disc.		0
Control 2	Autoclaved rice-complete disc(1:1)		0
<i>F. SCUTINELLUM</i> :			
01A	Grassfield	Lakeville	no
012C	Cultivated grass	Lakeville	no-0.50
022A	Cultivated grass	near Brownsville	no-0.50
022B	Cultivated grass	near Brownsville	no-0.50
022C	Cultivated grass	near Brownsville	no-0.50
022D	Cultivated grass	near Brownsville	no
022E	Cultivated grass	near Brownsville	no
022F	Cultivated grass	near Brownsville	no
022G	Cultivated grass	near Brownsville	no
025A	Grass	Farmstead	no-0.50
025B	Grass	Farmstead	no-0.50
025C	Grass	Farmstead	no-0.50
025D	Grass	Farmstead	no
025E	Grass	Farmstead	no
025F	Grass	Farmstead	no
025G	Grass	Farmstead	no
025H	Grass	Farmstead	no
025I	Grass	Farmstead	no
025J	Grass	Farmstead	no
025K	Grass	Farmstead	no
025L	Grass	Farmstead	no
025M	Grass	Farmstead	no
025N	Grass	Farmstead	no
025O	Grass	Farmstead	no
025P	Grass	Farmstead	no
025Q	Grass	Farmstead	no
025R	Grass	Farmstead	no
025S	Grass	Farmstead	no
025T	Grass	Farmstead	no
025U	Grass	Farmstead	no
025V	Grass	Farmstead	no
025W	Grass	Farmstead	no
025X	Grass	Farmstead	no
025Y	Grass	Farmstead	no
025Z	Grass	Farmstead	no

Table 2. (Continued)

<i>Fusarium</i> species and code number	Source	Area of origin	Mycelium*
<i>F. SCUTINELLUM</i> :			
031A	Grass	Oyster Bay	no
031B	Grass	Oyster Bay	0
031C	Grass	Alto	no
031D	Grass	Alto	no
031E	Grass	Alto	no
031F	Grass	Alto	no
031G	Grass	Alto	no
031H	Grass	Alto	no
031I	Grass	Alto	no
031J	Grass	Alto	no
031K	Grass	Alto	no
031L	Grass	Alto	no
031M	Grass	Alto	no
031N	Grass	Alto	no
031O	Grass	Alto	no
031P	Grass	Alto	no
031Q	Grass	Alto	no
031R	Grass	Alto	no
031S	Grass	Alto	no
031T	Grass	Alto	no
031U	Grass	Alto	no
031V	Grass	Alto	no
031W	Grass	Alto	no
031X	Grass	Alto	no
031Y	Grass	Alto	no
031Z	Grass	Alto	no

Table 2. (Continued)

<i>Fusarium</i> species and code number	Source	Area of origin	Mycelium*
<i>F. SCUTINELLUM</i> :			
031A	Grass	Alto	no-0.50
031B	Grass	Alto	no-0.50
031C	Grass	Alto	no-0.50
031D	Grass	Alto	no-0.50
031E	Grass	Alto	no-0.50
031F	Grass	Alto	no-0.50
031G	Grass	Alto	no-0.50
031H	Grass	Alto	no-0.50
031I	Grass	Alto	no-0.50
031J	Grass	Alto	no-0.50
031K	Grass	Alto	no-0.50
031L	Grass	Alto	no-0.50
031M	Grass	Alto	no-0.50
031N	Grass	Alto	no-0.50
031O	Grass	Alto	no-0.50
031P	Grass	Alto	no-0.50
031Q	Grass	Alto	no-0.50
031R	Grass	Alto	no-0.50
031S	Grass	Alto	no-0.50
031T	Grass	Alto	no-0.50
031U	Grass	Alto	no-0.50
031V	Grass	Alto	no-0.50
031W	Grass	Alto	no-0.50
031X	Grass	Alto	no-0.50
031Y	Grass	Alto	no-0.50
031Z	Grass	Alto	no-0.50

Table 2. (Continued)

<i>Fusarium</i> species and code number	Source	Area of origin	Mycelium*
<i>F. SCUTINELLUM</i> :			
031A	Cultivated grass area	Lakeville	no-0.50
031B	Cultivated grass area	Lakeville	no-0.50
031C	Cultivated grass area	Lakeville	no-0.50
031D	Cultivated grass area	Lakeville	no-0.50
031E	Cultivated grass area	Lakeville	no-0.50
031F	Cultivated grass area	Lakeville	no-0.50
031G	Cultivated grass area	Lakeville	no-0.50
031H	Cultivated grass area	Lakeville	no-0.50
031I	Cultivated grass area	Lakeville	no-0.50
031J	Cultivated grass area	Lakeville	no-0.50
031K	Cultivated grass area	Lakeville	no-0.50
031L	Cultivated grass area	Lakeville	no-0.50
031M	Cultivated grass area	Lakeville	no-0.50
031N	Cultivated grass area	Lakeville	no-0.50
031O	Cultivated grass area	Lakeville	no-0.50
031P	Cultivated grass area	Lakeville	no-0.50
031Q	Cultivated grass area	Lakeville	no-0.50
031R	Cultivated grass area	Lakeville	no-0.50
031S	Cultivated grass area	Lakeville	no-0.50
031T	Cultivated grass area	Lakeville	no-0.50
031U	Cultivated grass area	Lakeville	no-0.50
031V	Cultivated grass area	Lakeville	no-0.50
031W	Cultivated grass area	Lakeville	no-0.50
031X	Cultivated grass area	Lakeville	no-0.50
031Y	Cultivated grass area	Lakeville	no-0.50
031Z	Cultivated grass area	Lakeville	no-0.50

Table 2. (Continued)

<u>Fusarium species</u> and code number	Source	Area of origin	Mycoflora
<u>F. sporisorium:</u>			
063	Overwatered barley	O. Tocco	0
064A	Overwatered barley	O. Tocco	No
064B	Overwatered barley	O. Tocco	No
055A	Overwatered barley	O. Tocco	no-0, no
055B	Overwatered barley	O. Tocco	no-0, no
067B	Overwatered barley	same	0-1
000A	Everlasting	Oppland County	0-1
<u>F. subglutinans:</u>			
037A	Barley	Alta	0-1
037B	Barley	Alta	0-1
037C	Barley	Alta	0-1
037D	Barley	Alta	0-1
037E	Barley	Alta	0
037F	Barley	Alta	0
037G	Barley	Alta	0
037H	Barley	Alta	0
037I	Barley	Alta	0
037J	Barley	Alta	0
037K	Barley	Alta	0
037L	Barley	Alta	0
037M	Barley	Alta	0
037N	Barley	Alta	0
037O	Barley	Alta	0
037P	Barley	Alta	0
037Q	Barley	Alta	0
037R	Barley	Alta	0
037S	Barley	Alta	0
037T	Barley	Alta	0
037U	Barley	Alta	0
037V	Barley	Alta	0
037W	Barley	Alta	0
037X	Barley	Alta	0
037Y	Barley	Alta	0
037Z	Barley	Alta	0

Table 2. (Continued)

<u>Fusarium species</u> and code number	Source	Area of origin	Mycoflora
<u>F. sporisorium:</u>			
067A	Barley field	O. Tocco	0-1
067B	Barley field	O. Tocco	0-1
067C	Barley field	O. Tocco	0-1
067D	Barley field	O. Tocco	0-1
067E	Barley field	O. Tocco	0-1
067F	Barley field	O. Tocco	0-1
067G	Barley field	O. Tocco	0-1
067H	Barley field	O. Tocco	0-1
067I	Barley field	O. Tocco	0-1
067J	Barley field	O. Tocco	0-1
067K	Barley field	O. Tocco	0-1
067L	Barley field	O. Tocco	0-1
067M	Barley field	O. Tocco	0-1
067N	Barley field	O. Tocco	0-1
067O	Barley field	O. Tocco	0-1
067P	Barley field	O. Tocco	0-1
067Q	Barley field	O. Tocco	0-1
067R	Barley field	O. Tocco	0-1
067S	Barley field	O. Tocco	0-1
067T	Barley field	O. Tocco	0-1
067U	Barley field	O. Tocco	0-1
067V	Barley field	O. Tocco	0-1
067W	Barley field	O. Tocco	0-1
067X	Barley field	O. Tocco	0-1
067Y	Barley field	O. Tocco	0-1
067Z	Barley field	O. Tocco	0-1

Table 2. (Continued)

<u>Fusarium species</u> and code number	Source	Area of origin	Mycoflora
<u>F. sporisorium:</u> (Continued)			
067A	Barley field	O. Tocco	0-1
067B	Barley field	O. Tocco	0-1
067C	Barley field	O. Tocco	0-1
067D	Barley field	O. Tocco	0-1
067E	Barley field	O. Tocco	0-1
067F	Barley field	O. Tocco	0-1
067G	Barley field	O. Tocco	0-1
067H	Barley field	O. Tocco	0-1
067I	Barley field	O. Tocco	0-1
067J	Barley field	O. Tocco	0-1
067K	Barley field	O. Tocco	0-1
067L	Barley field	O. Tocco	0-1
067M	Barley field	O. Tocco	0-1
067N	Barley field	O. Tocco	0-1
067O	Barley field	O. Tocco	0-1
067P	Barley field	O. Tocco	0-1
067Q	Barley field	O. Tocco	0-1
067R	Barley field	O. Tocco	0-1
067S	Barley field	O. Tocco	0-1
067T	Barley field	O. Tocco	0-1
067U	Barley field	O. Tocco	0-1
067V	Barley field	O. Tocco	0-1
067W	Barley field	O. Tocco	0-1
067X	Barley field	O. Tocco	0-1
067Y	Barley field	O. Tocco	0-1
067Z	Barley field	O. Tocco	0-1

0 = no mycelia detected; 0-1 = unknown non-epitheliomycetes  
 mycelia produced by Fusarium spp. (M.A. 226, M.A. 181-182°C);  
 No = nonilluminated; F-c = Fusaria c; 0-1 = unknown non-  
 epitheliomycetes mycelia produced by Fusarium spp. (M.A. 400, M.A.  
 221-222°C); and F-0 = sporulation.

Table 3. Mycotoxins detected in extracts of fungal isolates used in these studies.

<u>Fusarium</u> species	Number isolates	Number of isolates producing mycotoxins					
		Trichothecenes	F-2	Fusarin c	YDP-1	HM-8*	M-1**
<u>F. acuminatum</u>	25	0	0	0	0	13	0
<u>F. avenaceum</u>	28	0	0	1	0	14	1
<u>F. culmorum</u>	01	0	1	0	0	0	0
<u>F. oxysporum</u>	31	0	0	0	0	14	12
<u>F. sambucinum</u>	40	0	0	0	0	2	34
Total no. of isolates	125	0	1	1	0	43	47

\* Metabolite of unknown structure produced by Fusarium spp. (M.W. 226, M.P. 181-182°C).

\*\* Metabolite of unknown structure produced by Fusarium spp. (M.W. 428, M.P. 222-223°C).

Code #	Fusarium Species	Area of origin	Death of rat within the feeding test	DON	3 ac- DON	Toxin		P-2	Monilic.	H-1
						H-8	P-X			
2 P <sub>63</sub> P <sub>9</sub>	<i>F. avenaceum</i>	Pasture leaves	1	x?			x		x	
1 P <sub>2</sub>	<i>F. culmorum</i>	South Auckland	1	x			x	x		
9 P <sub>10</sub>	<i>F. culmorum</i>	South Auckland	3	x			x	x		
6 P <sub>4</sub>	<i>F. oxysporum</i>	South Auckland	3	x	x			x		
8 P <sub>16</sub>	<i>F. oxysporum</i>	South Auckland	1	x		x <sup>a</sup>		x <sup>a</sup>	x	
9 P <sub>4</sub>	<i>F. oxysporum</i>	South Auckland	-					x	x	
9 P <sub>6</sub>	<i>F. sambucinum</i>	South Auckland	1	x			x	x		
3 A	<i>F. sambucinum</i>	Tara-high wet stocking	1							x
4 A	<i>F. sambucinum</i>	Tara-shade	1							x
5 A	<i>F. sambucinum</i>	Tara-shade	3	x?						x
6 A	<i>F. sambucinum</i>	Tara-random	1							x <sup>aa</sup>
5 B	<i>F. sambucinum</i>	Tara-shade	3	x?						x
1 D	<i>F. sambucinum</i>	Tara-shade	1							x
MS 1 A	<i>F. sambucinum</i>	Pukete Rd	1					x <sup>aa</sup>	x	
MS 4 A	<i>F. sambucinum</i>	Pukete Rd	-			x?	x	x <sup>a</sup>		

TABLE 4. Mycotoxins found in fifteen toxic isolates of *Fusarium* collected in New Zealand.



The extraction of the culture material was done with MeOH/H<sub>2</sub>O, ratio 55:45.  
The toxins were only confirmed by DC (TLC plates - Kieselgel 60, Merck 20 x 20 cm layer thickness 0.25 mm).

Developing systems were - chloroform/Methanol, ratio 9:1 and 4:1  
- chloroform/Aceton, ratio 3:1

Spraying reagent were - p - anisaldehyd  
- 20% H<sub>2</sub>SO<sub>4</sub>  
- 4 - (4 - Nitrobenzyl) purdine together with tetraethylene - pentamino.  
- 20% AL Cl<sub>3</sub> · H<sub>2</sub>O  
- 0.32% (w/v) aqueous 2,4 dinitrophenylhydrazine

Standards used for comparison: MAS, DAS, T-2, acethyl - T-2, HT-2, T-2 tetraol, DON, 3 acetyl-DON, 15 acetyl - DON, P-H, Nivalenol, HM-8, acetyl - HM-8, P-2, Moniliformin, H-1, TDP -1

Code:	x7	± not sure
	x*	± much
	x**	± very much

## DISTRIBUTION LIST

5 copies	Commander US Army Medical Research Institute of Infectious Diseases ATTN: SGRD-UIZ-M Fort Detrick, Frederick, MD 21701-5011
1 copy	Commander US Army Medical Research and Development Command ATTN: SGRD-RMI-S Fort Detrick, Frederick, MD 21701-5012
12 copies	Defense Technical Information Center (DTIC) ATTN: DTIC-DDAC Cameron Station Alexandria, VA 22304-6145
1 copy	Dean School of Medicine Uniformed Services University of the Health Sciences 4301 Jones Bridge Road Bethesda, MD 20814-4799
1 copy	Commandant Academy of Health Sciences, US Army ATTN: AHS-CDM Fort Sam Houston, TX 78234-6100